

SHIVAM INSTRUMENTS

(An ISO 9001:2008 certified company)



LABORATORY PULVERIZER(Hammer Mill Type)

This type of pulverizer is hammer Type fitted with a high speed motor. For fine powder of 200 mesh. It is fabricated on mild steel base fitted with forged steel hammers and complete with hopper, product container with lid & cloth product delivery tube.

TECHNICAL DATA:-

S.NO	DESCRIPTION	COMPLIANCE
1	MOC of outer body	Mild stel frame painted
2	Moc of hammer,liner body	Megnese steel /Stainless steel
3	Dia of drum	10"
4	Motor	3HP 440 volts 3 phase Perfect make
5	Starter	CG/Bentex/equivlent push button type starter
6	Input size	6-12.5mm Max.dry lime stone
7	Out put size	212 micron/72 mesh
8	Capacity	10-15KG /Hr. depending on the quality of coal.
9	Material Retention	2-5% Of feed material