

TEST SIEVES

Testing Sieves are perfectly suitable for laboratory and industrial uses where high quality sieving analysis is required. These test sieves are available in various standard sizes and can be customized to meet special requirements.

•TEST TIEVES OF COURSE SERIES:-

Made of perforated plate sheet with square holes. Available in 200 /300/450mm, dia. Aperture size:- 125,160,100,90,80,75,63,53,50,45,40,37.5,31.5,25,22.4,19,16,13.2,12.2,11.2, 10.2,10,95,8,6,7,6.3, 5.6,4.75,2.8mm

TEST SIEVES FINE SERIES IN 200MM DIA BRASS FRAME

BSS	ASTM	ISI	MICRONS	MM
4	5	--		4.75
5	6	340	1812	3.35
6	7	280	2411	2.80
7	8	240	2207	2.36
8	10	200	2057	2.00
10	12	170	1698	1.70
12	14	140	1005	1.40
14	16	120	1200	1.18
16	18	100	1000	1.00
18	20	85	850	0.850
22	25	70	710	0.710
25	30	60	600	0.600
30	35	50	500	0.500
36	40	40	425	0.425
44	45	35	355	0.355
52	50	30	300	0.300
60	60	25	250	0.250
72	70	20	212	0.212
85	80	18	180	0.180
100	100	15	150	0.150
120	120	12	125	0.125
150	140	10	106	0.106
170	170	9	90	0.90
200	200	8	75	0.75
240	230	7	63	0.63
300	270	6	53	0.53
325	300	5	45	0.45
350	325	4	44	0.044
400	380	3	40	0.040
500	---	--	25	0.025



TEST SIEVES



SHIVAM INSTRUMENTS
DELHI-110052
9312838091, 9811362676



mverma746@gmail.com | info@shivamlabinstruments.com
www.shivamlabinstruments.com